

**VWR**   
INTERNATIONAL

## Spectronic GENESYS 10 Bio UV-Visible Spectrophotometer



### Versatile, accurate UV-Vis for life science applications

In busy life-science labs, nothing can be left to chance. There's no room for error. Time is critical. Thermo Electron Corporation, with over 60 years of experience in UV-Visible spectrophotometry, has designed the Spectronic™ GENESYS™ 10 Bio spectrophotometer – the instrument you can count on to meet the demands of your laboratory.

#### Rapid, Reliable Nucleic Acid and Protein Analysis Results

The Spectronic GENESYS 10 Bio spectrophotometer puts leading-edge technology into a compact, convenient design. You get all the versatility and flexibility you need to obtain accurate, reliable measurements for your biotech research. The instrument features a user-friendly interface, powerful optics, rugged design, versatile software capabilities, optional built-in printer, and much more.

The time-saving software makes your job easier. Just put in your sample and obtain accurate results in seconds. Pre-programmed procedures enable even first-time users to:

- Measure DNA ratio and concentration with and without scanning
- Measure protein concentrations using Bradford, Lowry, BCA, and Biuret reagents
- Measure concentrations of dsDNA, ssDNA, RNA, and oligonucleotides
- Monitor cell growth
- Perform oligo calculations of molecular weight and T<sub>m</sub>

#### Plus, All the Capabilities of a Standard Lab Spectrophotometer

The Spectronic GENESYS 10 Bio spectrophotometer is not just dedicated to DNA and protein analyses. It is an ideal general-purpose spectrophotometer for such life science applications as multi-wavelength and simple kinetic-reaction measurements.

This instrument is monochromator-based, offering you a full wavelength range and superb accuracy.

#### Everything a Life Sciences Lab Needs and More...

- Customize your instrument for your specific applications. Use the unique SmartStart™ feature to easily create a power-up menu that displays only the tests that your lab uses. This saves time and reduces error.
- Save time by using the built-in Oligo calculator to determine the molecular weight and the theoretical T<sub>m</sub> (used to determine the PCR annealing temperature) of your oligonucleotides.
- Save time and reduce setup error by saving your programmed methods, and protect your saved methods with passwords.
- Easily ensure that your instrument is performing to specification with logs and documented evidence using built-in and/or accessory traceable instrument performance standards.

### Hardware

- Optical specifications comparable to instruments twice their size
- High-resolution graphics display
- Small footprint saves bench space
- Easy-to-use keyboard to set up tests, print results, and move cell changer

### Additional Features and Accessories

- Test tube holders
- Affordable sipper accessories
- Multi-cell holders measure multiple samples with the push of a button
- Micro-cells for small volumes
- Built-in printers show all parameters, date, time, test name, data, final results, and statistics

### For Busy Lab Operators, it is the Smart Choice

This spectrophotometer features the unique SmartStart screen that lets you select only the tests you require. Simply pre-program one, two, or more of your frequently used methods. When you turn the instrument on, the boot-up screen immediately shows the selected tests – without additional scrolling. After your data is acquired, you can print all instrument parameters, date, time, sample number, and data on an optional printer.


### DNA/Protein Concentration and DNA Purity

- Test results appear on the screen within seconds.
- Parameters for determining DNA purity, DNA concentration, and protein concentration are pre-selected; use default parameters or customize them to specific needs.
- Choose whether or not to display protein concentration on the results screen.
- An optional reference wavelength corrects for the effect of turbidity.
- Built-in dilution-factor calculator delivers direct concentration results.

DNA (260/280) 12:06pm 01Jul04 Cell #5			
ID#	Abs 260nm	Abs 280nm	Abs 320nm
1	0.725	0.410	0.010
	DNA Ratio =		1.79
	DNA Conc =		30.57 µg/mL
	Protein Conc =		79.33 µg/mL
2	0.294	0.162	0.005
	DNA Ratio =		1.84
	DNA Conc =		12.53 µg/mL
	Protein Conc =		24.80 µg/mL

### DNA/Protein Concentration and DNA Purity with Scan

- Provides an easy, one-step operation
- Displays scan of DNA sample and results
- Offers the option to add reference wavelength and dilution factor




### Direct UV/DNA Measurements

- Measures direct ssDNA, dsDNA, RNA, and oligo concentrations using absorbance at 260 nm
- Measures oligo concentrations using base sequence calculator
- Enables operators to select pre-set units or create new ones

Oligos (calc factor) 12:31pm 01Jul04 Cell #1			
ID#	Abs 260nm	Oligos µg/mL	Oligos pmol/µL
1	0.736	21.79	1.342
2	0.309	9.146	0.562
3	0.452	13.38	0.824

### Protein Concentration

- Offers seven standard curve methods for measuring protein – Bradford, Lowry, BCA (all with standard and micro) plus Biuret
- Allows use of pre-set or user-defined wavelength, units of measure, and curve fit
- Provides three UV methods – direct 280 nm, direct 205 nm, and Warburg-Christian calculation
- Offers print-outs that show the standard curve and results



Std Curve - Standards 12:06pm 01Jul04 Cell #5		
Std#	µg/mL	Abs 595nm
1	200	0.221
2	400	0.464
3	600	0.521
4	800	0.650
5	1000	0.883

Curve Fit = Linear Regr 2nd Order  
a0 = 0.00 a1 = 0.001  
a2 = 0.121  
Std Dev = 0.068


### Cell Growth

- Monitors cell growth by measuring absorbance at 600 nm
- Determines when to harvest cells

### Oligo Calculator

- Enter in a base sequence and the instrument calculates molecular weight, theoretical Tm, and oligo concentration factor.
- Choose from one of four Tm calculation methods.

Oligos 12:06pm 01Jul04	
Base Sequence =	ATCGTCGATTGAGCATCAGCATGACTAGATCAGATCGCG
Number of bases =	40
%GC =	47.56
DNA Mol. Wt. =	12388
DNA $\epsilon$ (260) =	441550
Conversion Factor =	28.06



Typical printout showing standard curve data and plot.

### Built to Perform. Built to Last

As a researcher, you need quick answers from your tests or data to support theories of molecular-reaction mechanisms. Whatever your need or application, the Spectronic GENESYS 10 Bio spectrophotometer can meet it.

The instrument features high-quality optics for accurate and precise measurements. Its long-life xenon lamp delivers an outstanding signal-to-noise ratio, negligible drift and no warm-up time for accurate, precise measurements over the 190 – 1100 nm wavelength range. It also provides calibrated wavelength peaks that are used as built-in wavelength standards for performance verification.

### Spectronic GENESYS 10 Bio Spectrophotometer – Additional Benefits

- 5 nm bandwidth for DNA peak resolution
- 100 nm range survey scan to find peak wavelengths
- Good stray-light rejection for excellent linearity
- Ability to save methods to non-volatile memory
- Customize your instrument using the SmartStart feature to display your tests on the boot-up menu
- Automatic cell correction for cell-to-cell precision with unmatched cuvettes

### VISION/ite Software for Expanded Capabilities

VISION/ite™ software offers a wide range of data handling, recording, and processing capabilities. You can easily perform wavelength scans, kinetics, single and multi-wavelength, and standard curve analyses with this software.

#### Scan – Collection of Sample Spectrum

This software application records a spectrum of your sample for measurement or characterization. Scan options include:

- A selectable wavelength range allowing full or specific spectrum investigations
- Peak pick which identifies wavelength peaks or valleys for further fixed-wavelength analysis
- Up to 13 spectral overlays
- Auto scaling
- Data saving
- Data file exporting

#### Rate – Measurements of Reaction Kinetics

This software application measures the change in your samples at a specific wavelength over a period of time. Rate capabilities include:

- Measurements of up to 100 samples
- Up to 13 graph overlays
- Automatic enzyme activity calculation
- Data rescale
- Automatic data storage
- Import and export of results

#### Fixed – Measurements at Pre-selected Wavelengths

This application performs measurements from one to 31 wavelengths which can be performed with a reference wavelength for automatic correction of a shifting background absorbance for DNA or turbid sample analyses. Results can be automatically stored and exported.

#### Quant – Determination of Sample Concentration

This software application measures sample concentrations based on either a standard curve or an entered calibration factor. Standard curve graphs and tabular data, along with calibration curve parameters (slope and intercept) are viewed on one screen. Optional parameters include:

- Four curve-fit options for non-linear standard curves
- Measurements of up to 20 standards
- The choice of 11 pre-programmed units
- Auto data save
- Data exports



## Faster Measurements Using Sample Holder of Your Choice

Before the software can work its magic, you need an instrument versatile enough to accommodate a wide range of optional sample holders. The Spectronic GENESYS 10 Bio spectrophotometer gives you an impressive range of choices for fast, efficient sample measurements.


### Longpath Cell Holders

- Used to measure low absorbance samples without the time and error of physically concentrating the sample
- Used to increase sensitivity of readings
- Hold rectangular or cylindrical cells
- Hold up to three 50 mm pathlength cells in 6-Position Cell Holder or one 100 mm pathlength cell in Single Cell Platform



### Flow-Through Accessory

- Great for fast sample throughput of batch runs
- Optical glass or quartz flowcell
- Selectable tubing and pumping times accommodate sample volumes from 2–15 mL
- Option to return sample to vessel



### Filter Holders

- Used to measure transmittance and cut-off wavelengths in glass and filters
- Both spring-loaded and adjustable holders are available



### Funnel Flowcell

- Fast sample throughput
- Ideal for quick, inexpensive bulk sample processing
- Used with external vacuum pump or water aspiration system
- Used on Single Cell Platform



### Test Tube Holder

- Ideal for COD tubes and glassware from reagent kits
- Holds test tubes up to 25 mm diameter and up to 102 mm tall



MODEL	VWR Cat. Number
Spectronic GENESYS 10 Bio, single cell	57971-352
Spectronic GENESYS 10 Bio, single cell holder, printer	57971-354
Spectronic GENESYS 10 Bio, 6-position cell changer	57971-356
Spectronic GENESYS 10 Bio, 6-position cell changer, printer	57971-358

## Spectronic GENESYS 10 Bio Spectrophotometer

Spectral Bandwidth	5 nm
Optical System	Split-beam, grating-based, dual detectors
Lamp Source: Typical Lifetime	Xenon: 5 years
Wavelength: Range	190 – 1100 nm
Accuracy	± 1.0 nm
Repeatability	± 0.5 nm
Graphic Display	320 x 240 pixel LCD, 3.8" x 2.8"
Photometric: Range	0.3 – 125%T; -0.1 – 3.0A; 0 – 9999 C
Readout	Absorbance, Transmittance, Concentration
Accuracy	0.5% of reading or 5 mA, whichever is greater, up to 2A
Noise	< 1 mA at 0A; < 2 mA at 2A, peak-to-peak at 340 nm
Drift	< 1 mA/hour after warm-up
Stray Light	< 0.1%T at 220 and 340 nm
Standard Interface	Bi-directional RS-232C
Standard Cell Holder	1-position or 6-position automatic cell holder
Keypad	Membrane keypad
Software	DNA ratio/concentration and protein concentration with or without scanning Direct oligo concentration at 260 nm Direct protein at 280 nm and 205 nm Protein standard curves Cell growth Oligo calculator: absorptivity, molecular weight, factor and theoretical Tm Absorbance, %T, Concentration Standard curve Absorbance ratio Absorbance difference Multiwavelength Kinetics Survey Scan Performance Validation
Test Storage	Up to 40 sets of test parameters
Languages	Software, printout, and Operator's Manual: English, French, German, Spanish, Italian (user-selectable)
Printer (optional)	40-column internal graphical
Power Requirements	Selected automatically; 100 – 240 volts
Dimensions	330 W x 410 D x 235 H mm; (13" x 16" x 9")
Weight	8.6 Kg (19 lb.)
Warranty	1 year

**VWR**  <sup>TM</sup>  
**INTERNATIONAL**

**vwr.com 1.800.932.5000**

**Order from Over 750,000 Products**

©2004 VWR International, Inc. All rights reserved. Printed in U.S.A. 0704 5M Lit. No. 20402