



VWR™ 4-in-1 EZ Tote™

**4-in-1 EZ Tote™**  
**and Integral**  
**Wash Bottles**



*Main photo and left inset: 73520-838.*

*Right inset: 73520-870, 872*

## VWR™ 4-in-1 EZ Tote™

- Versatile bottle and supplies carrier

If you're looking for a safe, easy and economical way to transport chemicals and lab items from store or stock rooms back to your lab, then choose the 4-in-1 EZ Tote. It carries 500mL, 1L or 4L bottles or may be used as an open tote to carry safety glasses, lab notebooks and a variety of lab items. Keep chemical bottles stored on the bench top or floor in the 4-in-1 EZ Tote – they are ready to go when you need them. Tote folds to store flat and is a snap to assemble. Patent pending.

| Cat. No.  | Description    | Case Qty. | Case Price |
|-----------|----------------|-----------|------------|
| 73520-838 | 4-in-1 EZ Tote | 4         | \$79.60    |

## VWR™ Integral Wash Bottles

- Easy to squeeze
- Translucent polyethylene

Wash bottles with integrally molded tube. Dispensing tip is provided with polyethylene closure cap and can be cut back to increase flow. Special design allows for easy grip and stability.

| Cat. No.  | Description     | Cap, mL | Case Qty. | Case Price |
|-----------|-----------------|---------|-----------|------------|
| 73520-870 | Wash Bottle, PE | 250     | 36        | \$90.00    |
| 73520-872 | Wash Bottle, PE | 500     | 24        | 73.20      |

**VWR:**  
Providing value at competitive prices



**vwr.com** 1.800.932.5000

Order from Over 750,000 Products

### Sales & Inventory Locations

#### Pacific Northwest

Anchorage, AK  
Boise, ID  
Brisbane, CA  
Seattle, WA  
Tualatin, OR

#### Southwest

Albuquerque, NM  
Aurora, CO  
Phoenix, AZ  
San Diego, CA  
San Dimas, CA

#### Midwest

Batavia, IL  
Indianapolis, IN  
Minneapolis, MN  
St. Louis, MO

#### Gulf

Houston, TX  
Lake Charles, LA  
Sugarland, TX

#### Northeast

Marlborough, MA  
Rochester, NY

#### Mid-Atlantic

Bridgeport, NJ  
Clarksburg, MD

#### Southeast

Morrisville, NC  
Suwanee, GA

#### Canada

Mississauga, Ontario  
Ville Mont-Royal,  
Québec

### Call Direct for Specialized Service Locations

#### International

911 Commerce Ct.  
Buffalo Grove, IL 60089  
847-520-6170  
fax: 847-229-0788

#### Puerto Rico

Carr. #869, Km. 1.5 M-4  
Royal Industrial Park  
Catano, PR 00962  
800-932-5000  
787-275-4500  
fax: 787-788-4320

#### Global and

#### Corporate Accounts

1310 Goshen Pkwy.  
West Chester, PA 19380  
610-431-1700  
fax: 610-431-9174

#### Mexico

Super Av. Lomas Verdes  
No. 464 2o Piso Col.  
Los Alamos C.P. 53230  
Naucalpan de Juarez,  
Edo, de Mexico 2122-1600  
011-52-552-122-1630  
fax: 011-52-555-344-7532

#### Production Supplies and Services

975 Overland Ct.  
San Dimas, CA 91773  
877-VWR-PSCS  
(877-897-7727)  
fax: 888-890-9124

#### Furniture Group

3021 Gateway Drive  
Suite 280  
Irving, TX 75063  
888-624-2432  
fax: 972-753-1389

#### VWR Laboratory Distributors Services

911 Commerce Ct.  
Buffalo Grove, IL 60089  
800-873-8971  
fax: 847-229-0668

Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.

©2003 VWR International, Inc. All rights reserved. Printed in U.S.A.