



vwr.com
1.800.932.5000

Order from Over 400,000 Products

Sales & Inventory Locations

Pacific Northwest
Anchorage, AK
Boise, ID
Brisbane, CA
Seattle, WA
Tualatin, OR

Southwest
Albuquerque, NM
Aurora, CO
Phoenix, AZ
San Diego, CA
San Dimas, CA

Midwest
Batavia, IL
Indianapolis, IN
Minneapolis, MN
St. Louis, MO

Gulf
Houston, TX
Lake Charles, LA
Sugarland, TX

Northeast
Marlborough, MA
Rochester, NY

Mid-Atlantic
Bridgeport, NJ
Clarksburg, MD

Southeast
Miami, FL
Morrisville, NC
Oak Ridge, TN
Suwanee, GA

Canada
Mississauga, Ontario
Ville Mont-Royal,
Québec

Call Direct for Specialized Service Locations

International
911 Commerce Ct.
Buffalo Grove, IL 60089
847-520-6170
fax: 847-229-0788

Puerto Rico
Carr. #869, Km. 1.5 M-4
Royal Industrial Park
Catano, PR 00962
800-932-5000
787-275-4500
fax: 787-788-4320

Global and Corporate Accounts
1310 Goshen Pkwy.
West Chester, PA 19380
610-431-1700
fax: 610-431-9174

Production Supplies and Services
975 Overland Ct.
San Dimas, CA 91773
877-VWR-PSCS
(877-897-7727)
fax: 888-890-9124

Furniture Division
3021 Gateway Drive
Suite 280
Irving, TX 75063
888-624-2432
fax: 972-753-1389

VWR Laboratory Distributors Services
911 Commerce Ct.
Buffalo Grove, IL 60089
800-873-8971
fax: 847-229-0668

Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.
©2009 VWR International, Inc. All rights reserved. Printed in U.S.A. 05K02



VWR™ Autoclavable Bags

**Choose VWR
Autoclavable
Bags for
Performance
and Safety**



VWR Catalog number

VRW 1.5 mil Autoclavable Bags



VWR autoclave bags provide an efficient, safe solution for the messy and potentially dangerous process of disposing of biohazardous waste material.

High strength, high temperature resistant polypropylene construction allows perfect

adaptability to standard laboratory disposal routines, including steam autoclaving. Rated for use at 123 C (275 F).

These bags are available clear or red with proper universal biohazard warnings and use instructions, or plain clear bags with no printing.

Made with Suprop brand PP by Tufpak, Inc.

	Dimensions , W x H cm(in)	Dimensions , W x H in	Catalog No.	Case Of/Cost
Red Biohazard Bags	20.3 x 30.5	(8 x 12)	14220-000	1000/195.65 100/pk 19.57
	30.5 x 61	(12 x 24)	14220-002	1000/389.45 100/pk 38.95
	48.3 x 61	(19 x 24)	14220-004	500/374.45 100/pk 74.89
	61 x 91.4	(24 x 36)	14220-006	500/611.00 100/pk 122.20
Clear Biohazard Bags	25.4 x 38.1	(10 x 15)	14220-008	1000/228.70
	30.5 x 61	(12 x 24)	14220-010	1000/349.23
	61 x 76.2	(24 x 30)	14220-012	500/423.07
	61 x 91.4	(24 x 36)	14220-014	500/499.24
Plain Clear Autoclave Bags	20.3 x 30.5	(8 x 12)	14220-016	1000/129.30 100/pk 12.93
	25.4 x 38.1	(10 x 15)	14220-018	1000/149.30 100/pk 14.93
	30.5 x 61	(12 x 24)	14220-020	1000/282.14
	61 x 76.2	(24 x 30)	14220-022	500/345.37
	61 x 91.4	(24 x 36)	14220-024	500/387.60

VRW 2.0 mil Autoclavable BIOHAZARD Bags



Ideal for processing of biohazardous materials which require autoclaving prior to disposal. Available with and without temperature indicator patch. High integrity polypropylene bags withstand 135 C (275 F) maximum autoclave temperature. Twist-tie

closure allows bag to expand and vent steam during sterilization. Bags design resist puncture, tears, and leaks during proper handling. Meets the 165 gram dart drop standard as defined by ASTM D1709-98.

Available in red, orange or clear bags. Each imprinted with the universal biohazard symbol.

Bags with temperature indicator patch have wording " CONTENTS HAVE BEEN AUTOCLAVED " printed with a chemically active ink. Before autoclaving the ink is clear on a white background, after auto-claving words appear in a dark brown color. Providing proof that bag has been exposed to high temperatures for sterilization.

Made with Suprop brand PP by Tufpak, Inc.

	Dimensions , W x H	cm (in.)	Catalog No.	Case Of/Cost
Clear Biohazard Bags	30.5 x 61	(12 x 24)	14220-066	200/112.21
	61 x 76.2	(24 x 30)	14220-068	200/180.04
	61 x 91.4	(24 x 36)	14220-070	200/187.82
Red Biohazard Bags	20.3 x 30.5	(8 x 12)	14220-056	400/90.45
	35.6 x 48.3	(14 x 19)	14220-058	200/91.43
	48.3 x 61	(19 x 24)	14220-060	200/153.05
	63.5 x 88.9	(25 x 35)	14220-062	200/244.11
	91.4 x 114.3	(36 x 45)	14220-064	100/226.62
Red Biohazard Bags w/Temp. Indicator	20.3 x 30.5	(8 x 12)	14220-026	400/127.51
	35.6 x 48.3	(14 x 19)	14220-028	200/117.38
	48.3 x 61	(19 x 24)	14220-030	200/159.95
	63.5 x 88.9	(25 x 35)	14220-032	200/263.59
	91.4 x 114.3	(36 x 45)	14220-034	100/268.26
Orange Biohazard Bags w/Temp. Indicator	20.3 x 30.5	(8 x 12)	14220-072	400/128.48
	35.6 x 48.3	(14 x 19)	14220-074	200/84.34
	48.3 x 61	(19 x 24)	14220-076	200/152.53
	63.5 x 88.9	(25 x 35)	14220-078	200/232.43
	91.4 x 114.3	(36 x 45)	14220-080	100/195.04
Clear Autoclave Bags	20.3 x 30.5	(8 x 12)	14220-046	400/149.66
	35.6 x 48.3	(14 x 19)	14220-048	200/114.85
	48.3 x 61	(19 x 24)	14220-050	200/155.85
	63.5 x 88.9	(25 x 35)	14220-052	200/243.67
	91.4 x 114.3	(36 x 45)	14220-054	100/270.31

VRW 2.0 mil Autoclavable Bags NONhazardous Waste



Intended for use in autoclaving nonhazardous materials and waste. Bags are printed with usage instructions and right-on block. Polypropylene bags withstand 135 C (275 F)

maximum autoclave temperature. Twist-tie closure allows bag to expand and vent steam during sterilization. Bags design resist puncture, tears, and leaks during proper handling. Meets the 165 gram dart drop standard as defined by ASTM D1709-98. Available in clear.

VRW 2.0 mil HIGH STRENGTH Autoclavable Bags



Engineered to meet the most demanding regulatory specifications known. They exceed the 480 gram tear test as defined by ASTM D1299-94A in both directions.

These biohazard bags meet all State requirements for hazardous waste handling. Including Florida and Pennsylvania Department Of Health Codes. Biohazard bags have temperature indicator.

Made with Suprop + brand PP by Tufpak, Inc.

	Dimensions , W x H	cm (in.)	Catalog No.	Case Of/Cost
Clear Autoclave Bags	20.3 x 30.5	(8 x 12)	14220-036	400/127.00
	35.6 x 48.3	(14 x 19)	14220-038	200/113.61
	48.3 x 61	(19 x 24)	14220-040	200/150.72
	63.5 x 88.9	(25 x 35)	14220-042	200/234.92
	91.4 x 114.3	(36 x 45)	14220-044	100/224.07

	Dimensions , W x H	cm (in.)	Catalog No.	Case Of/Cost
Red Biohazard Bags w/Temp. Indicator	20.3 x 30.5	(8 x 12)	NEW	400/
	30.5 x 61	(12 x 24)	14220-090	200/125.12
	48.3 x 61	(19 x 24)	NEW	200/
	61 x 76.2	(24 x 30)	14220-092	200/243.84
	61 x 91.4	(24 x 36)	14220-094	200/269.69
	63.5 x 88.9	(25 x 35)	14220-096	200/280.60
	96.5 x 121.9	(38 x 48)	14220-098	100/190.98
Red Autoclave Bags NO PRINT	30.5 x 61	(12 x 24)	14220-100	200/99.23
	61 x 76.2	(24 x 30)	14220-102	200/217.57
	61 x 91.4	(24 x 36)	14220-104	200/229.60

VRW 3.0 mil Double Thick Autoclavable Bags



High strength, high temperature resistant 3 mil polypropylene bags. Autoclavable up to 135 C (275 F). A heavy duty bag with extra puncture resistance and impact strength.

	Dimensions , W x H	cm (in.)	Catalog No.	Case Of/Cost
Red Biohazard Bags w/Temp. Indicator	30.5 x 61	(12 x 24)	14220-082	250/133.73
	61 x 76.2	(24 x 30)	14220-084	250/285.58
	61 x 91.4	(24 x 36)	14220-086	250/331.98
	96.5 x 122	(38 x 48)	14220-088	50/129.70

VRW 4.0 mil Double Thick Autoclavable Bags

4 mil high molecular high density polyethylene (HDPE) bags. Designed for autoclaving animal cages and other sharp edged items. Maximum autoclave temperature is 121 C (250 F). Ideal for non waste items.

Made with SecuraT brand PE by Tufpak, Inc.

	Dimensions , W x H	cm (in.)	Catalog No.	Case Of/Cost
Red Biohazard Bags w/Temp. Indicator	30.5 x 61	(12 x 24)	11215-822	250/125.70 50/pk 25.14
	61 x 76.2	(24 x 30)	11215-823	250/256.40 50/pk 51.28
	61 x 91.4	(24 x 36)	11215-826	250/388.65 50/pk 77.73
	96.5 x 122	(38 x 48)	11215-824	50/165.79