



VWR™ Beakers
VWR™ Specimen Containers
VWR™ Centrifuge Bottles
VWR™ Disposable Tissue Grinders

**Reaching the
Highest Levels
of Quality**





VWR™ Beakers

High-quality VWR beakers are the premier choice for scientific applications. All VWR beakers are transparent for clear content visibility. They have molded-in graduations for easily estimating volumes, as well as molded-in pouring spouts and rolled lip to prevent drips. VWR beakers are made of the highest quality polypropylene and are resistant to most commonly-used laboratory reagents, including a variety of strong acids, bases and organic solvents.

Cat. No.	Description	Graduation Range	Case Price	Pack Price
47751-690	15 mL Graduated Beaker	5–15 mL	1000/\$123.60	100/\$15.67
47751-692	30 mL Graduated Beaker	5–30 mL	1000/150.38	100/19.06
47751-694	50 mL Graduated Beaker	5–50 mL	1000/164.80	100/20.88
47751-696	100 mL Graduated Beaker	5–100 mL	1000/220.00	100/27.88
47751-698	150 mL Graduated Beaker	5–150 mL	1000/377.00	100/47.77
47751-700	250 mL Graduated Beaker	10–250 mL	500/229.00	50/29.02
47751-702	400 mL Graduated Beaker	25–400 mL	500/314.00	50/39.79
47751-705	600 mL Graduated Beaker	25–600 mL	500/361.00	25/22.86



VWR™ Specimen Containers

These lightweight, all-purpose specimen containers can be used for almost any liquid or solid specimen.

Molded of high-density polyethylene, (except for 32 oz., which is molded of polypropylene), containers are virtually indestructible in normal use.

Inert to most chemicals, including weak acids and bases, and resistant to freezing or boiling. The low, wide-mouth design provides ample access for filling container cleanly, with fewer spills, and makes it less likely to tip over. The tight-fitting cover keeps odors inside and minimizes risk of leakage.

Translucent walls keep contents visible. Economically priced to be disposable, yet reusable, too.

Cat. No.	Description	Case Price
44333-000	4 oz. Translucent Multipurpose Container	300/\$80.87
44333-002	8 oz. Translucent Multipurpose Container	100/36.88
44333-004	16 oz. Translucent Multipurpose Container	100/49.40
44333-006	32 oz. Translucent Multipurpose Container	100/79.44
44333-008	64 oz. Translucent Multipurpose Container	50/59.13
44333-010	86 oz. Translucent Multipurpose Container	25/38.26
44333-012	172 oz. Translucent Multipurpose Container	10/31.31



VWR™ Centrifuge Bottles with Caps, Spherical-Bottom

- Spherical-bottom bottles available in polycarbonate or polypropylene
- Resistant to most chemicals
- Wide-mouth design makes filling and decanting easy
- Choice of polypropylene screw cap or Noryl® sealing caps for a secure seal and easy gripping
- Autoclavable to 121°C at 15 psi for 20 minutes
- Large writing surface
- Printed volume graduations
- Precision molded construction

Polycarbonate Bottles

Cat. No.	Description	Max. RCF, g	Case Price
47735-692	250 mL Centrifuge Bottle w/screw cap	8250	36/\$275.83
47735-700	500 mL Centrifuge Bottle w/screw cap	8250	24/312.97
47735-694	250 mL Centrifuge Bottle w/seal cap	31,700	36/456.19
47735-702	500 mL Centrifuge Bottle w/seal cap	13,850	24/371.32

Polypropylene Bottles

Cat. No.	Description	Max. RCF, g	Case Price
47735-688	250 mL Centrifuge Bottle w/screw cap	8250	36/\$190.96
47735-696	500 mL Centrifuge Bottle w/screw cap	8250	24/153.83
47735-704	1000 mL Centrifuge Bottle w/screw cap	7400	16/196.27
47735-690	250 mL Centrifuge Bottle w/seal cap	27,500	36/386.25
47735-698	500 mL Centrifuge Bottle w/seal cap	13,500	24/350.10

Centrifuge Bottle Compatibility Guide

High-Speed Rotors	Beckman Rotors	Sorvall Rotors	Thermo IEC Rotors	VWR Cat. Number
6 x 500 mL Fixed Angle	JA-10	GS-3	878	47735-698; 47735-702
	JLA-10.500	SLA-3000	—	47735-690**; 47735-694**
6 x 250 mL Fixed Angle	JA-14	GSA*	877	47735-690; 47735-694
	JLA-16.250	SLA-1500	872	
4 x 250 mL Swinging Bucket	JS-7.5	HS-4	963	47735-690; 47735-694
	JS-4.3**	SH-3000**	947	

Low-Speed Rotors	Beckman Rotors	Sorvall Rotors	Thermo IEC Rotors	Juan Rotors	VWR Cat. Number
6 x 1000 mL Swinging Bucket	JS-4.2	H-6000A	966	C60	47735-704; 47735-698***
	JS-3.0	—	981	P60	47735-700***; 47735-688
		—	949	—	47735-692**
4 x 1000 mL Swinging Bucket 688**	JS-4.0	H-4000	256	P40	47735-704; 47735-696***
	JS-5.2	—	276	—	47735-700***; 47735-688***
					47735-692**
8 x 250 mL Swinging Bucket	—	H-2000B	—	—	47735-688; 47735-692
6 x 250 mL Swinging Bucket	—	—	259	—	47735-688; 47735-692
4 x 250 mL Swinging Bucket	GH-3.7**	H-6000B	219**; 218	M4	47735-688; 47735-692
	TH-4**	RTH-750**	210**; 284	T4	
			224; 268		
			245; 243		
8 x 250 mL Fixed Angle	—	—	850S	—	47735-688; 47735-692

* Requires Nalgene, Adapter 21020-201. **With rotor manufacturer's 250 mL adapters. ***With rotor manufacturer's 500 mL adapters.

The Centrifuge Bottle Compatibility chart is a guide to help you select the right VWR Brand Centrifuge bottle for your application and equipment. The information is believed to be reliable and accurate, but since conditions of use are beyond its control, VWR cannot guarantee rated performance in all centrifuges and rotors.

VWR™ Disposable Tissue Grinders

This tissue grinding system was designed with safety and ease-of-use in mind. The patented closed system has the ability to grind the specimen in a sealed container, helping to prevent any risk of contact with the sample. The molded-in abrasive surface does not leave any sediment to obstruct the sample when examining it. This system comes fully-assembled and Gamma-sterilized. An extra cap and patient identification label are provided for storage, transport, and added safety. The 35 mL grinder is designed to be free standing to allow for easier handling.



Cat. No.	Description	Case Price
47732-446	VWR 15 mL Conical Tissue Grinder–Sterile	10/\$86.47
47732-448	VWR 35 mL Free Standing Tissue Grinder–Sterile	10/86.47
47732-450	VWR 50 mL Conical Tissue Grinder–Sterile	10/86.47



vwr.com 1.800.932.5000

Order from Over 750,000 Products

Sales & Inventory Locations

Pacific Northwest

Anchorage, AK
Boise, ID
Brisbane, CA
Seattle, WA
Tualatin, OR

Southwest

Albuquerque, NM
Aurora, CO
Phoenix, AZ
San Diego, CA
San Dimas, CA

Midwest

Batavia, IL
Indianapolis, IN
Minneapolis, MN
St. Louis, MO

Gulf

Lake Charles, LA
Sugarland, TX

Northeast

Franklin, MA
Rochester, NY

Mid-Atlantic

Bridgeport, NJ
Clarksburg, MD

Southeast

Morrisville, NC
Suwanee, GA

Canada

Mississauga, Ontario
Ville Mont-Royal,
Québec

Call Direct for Specialized Service Locations

International

911 Commerce Ct.
Buffalo Grove, IL 60089
847-520-6170
fax: 847-229-0788

Puerto Rico

Carr. #869, Km. 1.5 M-4
Royal Industrial Park
Catano, PR 00962
800-932-5000
787-275-4500
fax: 787-788-4320

Global and Corporate Accounts

1310 Goshen Pkwy.
West Chester, PA 19380
610-431-1700
fax: 610-431-9174

Mexico

Super Av. Lomas Verdes
No. 464 2o Piso Col.
Los Alamos C.P. 53230
Naucalpan de Juarez,
Edo, de Mexico 2122-1600
011-52-552-122-1630
fax: 011-52-555-344-7532

Production Supplies and Services

975 Overland Ct.
San Dimas, CA 91773
877-VWR-PSCS
(877-897-7727)
fax: 888-890-9124

Furniture Group

3021 Gateway Drive
Suite 280
Irving, TX 75063
888-624-2432
fax: 972-753-1389

VWR Laboratory Distributors Services

911 Commerce Ct.
Buffalo Grove, IL 60089
800-873-8971
fax: 847-229-0668

VWR:

Providing value at competitive prices



Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.

©2004 VWR International, Inc. All rights reserved. Printed in U.S.A.

1/04

0404R (5M) Lit. No. 99869