



VWR Signature™ Carboys

VWR Signature
Carboys
Contain and
Deliver



73520-548, 73520-550, 73520-554



VWR Signature™ Carboys

- Carboys with and without spigots
- Distinctive hexagonal shape for extra strength
- Molded handgrips

VWR Signature Carboys offer all the best features for laboratory applications. VWR Signature Carboys come in three popular sizes: 5 L, 10 L and 20 L, with or without spigots. These HDPE containers have excellent chemical resistance and comply with FDA requirements for food use. Distinctive hexagonal shape provides extra strength and rigidity. Molded handgrips allow easier lifting and carrying. Wide 100 mm neck provides easy access for cleaning or filling. Standard carboys graduated to contain; carboys with spigots graduated to deliver.

Product Description

Cat. No.	Description	Case Qty.	Each	Case
73520-546	5L VWR Signature Carboy	6	\$33.84	\$203.04
73520-548	5L VWR Signature Carboy with spigot	6	\$49.36	\$296.16
73520-550	10L VWR Signature Carboy	6	\$47.84	\$287.04
73520-552	10L VWR Signature Carboy with spigot	6	\$63.36	\$380.16
73520-554	20L VWR Signature Carboy	4	\$60.32	\$241.28
73520-556	20L VWR Signature Carboy with spigot	4	\$75.84	\$303.36
73520-560	Carboy Replacement Spigot, 12 mm	12	\$7.50	\$90.00
73520-562	Carboy Replacement Closure, 100 mm	12	\$2.16	\$26.02



vwr.com 1.800.932.5000

Order from Over 400,000 Products

Sales & Inventory Locations

Pacific Northwest

Anchorage, AK
Boise, ID
Brisbane, CA
Seattle, WA
Tualatin, OR

Southwest

Albuquerque, NM
Aurora, CO
Phoenix, AZ
San Diego, CA
San Dimas, CA

Midwest

Batavia, IL
Indianapolis, IN
Minneapolis, MN
St. Louis, MO

Gulf

Houston, TX
Lake Charles, LA
Sugarland, TX

Northeast

Marlborough, MA
Rochester, NY

Mid-Atlantic

Bridgeport, NJ
Clarksburg, MD

Southeast

Miami, FL
Morrisville, NC
Oak Ridge, TN
Suwanee, GA

Canada

Mississauga, Ontario
Ville Mont-Royal,
Québec

Call Direct for Specialized Service Locations

International

911 Commerce Ct.
Buffalo Grove, IL 60089
847-520-6170
fax: 847-229-0788

Puerto Rico

Carr. #869, Km. 1.5 M-4
Royal Industrial Park
Catano, PR 00962
800-932-5000
787-275-4500
fax: 787-788-4320

Global and Corporate Accounts

1310 Goshen Pkwy.
West Chester, PA 19380
610-431-1700
fax: 610-431-9174

Production Supplies and Services

975 Overland Ct.
San Dimas, CA 91773
877-VWR-PSCS
(877-897-7727)
fax: 888-890-9124

Furniture Division

3021 Gateway Drive
Suite 280
Irving, TX 75063
888-624-2432
fax: 972-753-1389

VWR Laboratory Distributors Services

911 Commerce Ct.
Buffalo Grove, IL 60089
800-873-8971
fax: 847-229-0668

VWR Signature Products:

Delivering high quality and superior performance



Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.

©2003 VWR International, Inc. All rights reserved. Printed in U.S.A.

05/02

0103R (20M) Lit. No. 99907