



VWR™ SuperClear™ Centrifuge Tubes

For High
Performance
Centrifugation
and Storage
to 15,000* x g



Cat.#21008-103, 21008-178

VWR™ SuperClear™ Centrifuge Tubes for Worry-Free Centrifugation.

VWR™ SuperClear centrifuge tubes are available in 15mL or 50mL sizes that fit popular centrifuges. All printed tubes feature extra large solvent resistant white labeling areas and black graduations. Available in sterile or non sterile bags, clean patented collapsible racks, or in bulk packaging. Every sterile lot is tested and is free of endotoxin and nuclease.

Ordering Information



These tubes are so clear you can see the VWR logo in reverse on the opposite side of the tube.

15mL SuperClear™ Tubes, Spin to 15,000 x g*

Cat. No.	Description	Package	Price
21008-214	15mL Tubes	50/Bag, 10 bags/Case	500/\$150.38
21008-216	15mL Tubes	50/Bag, 10bags/Case, Sterile	500/164.34
21008-089	15mL Tubes	25/Rack, 2 Racks/Pkg 10 Pkg/Case	500/153.14
21008-103	15mL Tubes	25/Rack, 2 Racks/Pkg 10 Pkg/Case	500/170.95
21008-105	15mL Tubes	Bulk Packed, 500/Bag	500/142.06



These tubes are so clear you can see the VWR logo in reverse on the opposite side of the tube.

50mL SuperClear™ Tubes, Spin to 15,000 x g*

Cat. No.	Description	Package	Price
21008-240	50mL Tubes	50/Bag, 10 bags/Case	500/\$143.22
21008-242	50mL Tubes	50/Bag, 10bags/Case, Sterile	500/207.14
21008-169	50mL Tubes	25/Rack, 20 Pkg/Case	500/189.38
21008-178	50mL Tubes	25/Rack, 20 Pkg/Case, Sterile	500/220.33
21008-177	50mL Tubes	Bulk Packed, 500/Bag	500/159.14

*Organic solvents and other strong chemicals may reduce maximum R.C.F. Always use appropriate cushions and adapters.

VWR:
Providing value at competitive prices



VWR™ SuperClear™ tubes withstand breakage better than other major tube brands. The photo above shows the results of a test of market samples of tubes compared to VWR SuperClear tubes. All tubes were filled with 45mL of distilled water and spun unsupported for 15 minutes at 20,000 x g. All the other brands of tubes broke. Only VWR SuperClear tubes with the purple cap protected the sample and the centrifuge.

VWR SuperClear tubes are real centrifuge tubes designed for spinning not simply for storage. They can be spun reliably to 15,000 x g. Protect your samples and your work with VWR SuperClear tubes.

VWR 
INTERNATIONAL

vwr.com 1.800.932.5000

Order from Over 750,000 Products

Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.

©2004 VWR International, Inc. All rights reserved. Printed in U.S.A.

SuperClear™ is a trademark of Labcon North America. VWR and the VWR logo and design are trademarks of VWR International, Inc. The purple color cap for centrifuge tubes is a registered trademark of Labcon North America.

0504 10M Lit. No. 92011REV