




## VWR™ Colorfrost® Plus Slides

**Top-Quality**

**Colorfrost®**

**Plus Slides**



*From left to right: 48312-908, 48312-906, 48312-902, 48312-904, 48311-703*

## VWR™ Colorfrost® Plus Slides


- Saves Time & Money – Studies show that purchasing adhesion slides saves time and money over coating slides in the lab.
- Electrostatically attracts fresh frozen tissue sections and cytology preparations.
- No need to handle dangerous reagents used in treating slides in the lab.
- Eliminates background staining issues that may occur when using gelatin.
- Slides can be personalized for added identification providing permanence, convenience and cost savings.
- Permanent Colorfrost® tab stands up to alcohol and xylene.
- Plus slides have a longer shelf life than Poly-L-lysine.
- Superfrost® (white slides) and Colorfrost® (colored slides) are the highest quality color-coded microslides available.

Cat. No.	Description*	Gross	Case of Gross
48311-703	White, Plus, 25 x 75 - 1mm	\$72.25	\$614.10
48312-902	Blue, Plus, 25 x 75 - 1mm	65.16	553.89
48312-904	Green, Plus, 25 x 75 - 1mm	62.65	523.53
48312-906	Pink, Plus, 25 x 75 - 1mm	65.16	553.89
48312-908	Yellow, Plus, 25 x 75 - 1m	62.65	532.52

\*Other colors available upon request.

VWR:

Providing value at competitive prices



**vwr.com** 1.800.932.5000

**Order from Over 400,000 Products**

### Sales & Inventory Locations

#### Pacific Northwest

Anchorage, AK  
Boise, ID  
Brisbane, CA  
Seattle, WA  
Tualatin, OR

#### Southwest

Albuquerque, NM  
Aurora, CO  
Phoenix, AZ  
San Diego, CA  
San Dimas, CA

#### Midwest

Batavia, IL  
Indianapolis, IN  
Minneapolis, MN  
St. Louis, MO

#### Gulf

Houston, TX  
Lake Charles, LA  
Sugarland, TX

#### Northeast

Marlborough, MA  
Rochester, NY

#### Mid-Atlantic

Bridgeport, NJ  
Clarksburg, MD

#### Southeast

Miami, FL  
Morrisville, NC  
Oak Ridge, TN  
Suwanee, GA

#### Canada

Mississauga, Ontario  
Ville Mont-Royal,  
Québec

### Call Direct for Specialized Service Locations

#### International

911 Commerce Ct.  
Buffalo Grove, IL 60089  
847-520-6170  
fax: 847-229-0788

#### Puerto Rico

Carr. #869, Km. 1.5 M-4  
Royal Industrial Park  
Catano, PR 00962  
800-932-5000  
787-275-4500  
fax: 787-788-4320

#### Global and Corporate Accounts

1310 Goshen Pkwy.  
West Chester, PA 19380  
610-431-1700  
fax: 610-431-9174

#### Production Supplies and Services

975 Overland Ct.  
San Dimas, CA 91773  
877-VWR-PSCS  
(877-897-7727)  
fax: 888-890-9124

#### Furniture Division

3021 Gateway Drive  
Suite 280  
Irving, TX 75063  
888-624-2432  
fax: 972-753-1389

#### VWR Laboratory Distributors Services

911 Commerce Ct.  
Buffalo Grove, IL 60089  
800-873-8971  
fax: 847-229-0668

Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.

Colorfrost is a registered trademark of Erie Scientific Company. US Patent No. 4,624,882

©2002 VWR International, Inc. All rights reserved. Printed in U.S.A.

1002 (10M) Lit. No. 99888