



VWR™ CryoPro® Inventory Storage Boxes and Grids

Inventory Control
For Mechanical
and Cryogenic
Freezers



CryoPro



VWR™ CryoPro® Cryogenic Inventory Storage Boxes and Grids

LN₂ Cryogenic Freezer Boxes with drain slots

VWR Cat#	Style	Dimensions (In.)	List Pricing
82007-158	100 cell	5" x 5" x 2"	\$4.00/ ea.
82007-160	25 cell	2 3/4" x 2 3/4" x 2"	\$3.45/ ea.

LN₂ cryogenic freezer boxes are constructed of durable fiberboard with a protective moisture repellent coating. Boxes feature slots enabling safe liquid nitrogen drainage and reducing LN₂ consumption. Boxes come standard with dividers installed. Lids are printed with a numbered grid for superior inventory control.



82007-158, 82007-160

Mechanical Cryogenic Freezer Boxes (No drain slots - NOT FOR USE WITH LN₂)

VWR Cat#	Style	Dimensions (In.)	List Pricing
82007-142	no divider	5" x 5" x 2"	\$2.15/ ea.
82007-144	no divider	5" x 5" x 3"	\$2.80/ ea.
82007-162	81 cell divider	5" x 5" x 2"	\$2.35/ ea.

Mechanical Cryogenic Freezer Boxes are constructed of durable fiberboard with a protective moisture repellent coating. Box #82007-162, the most popular size, comes standard with an 81 Cell divider installed and numbered grid printed on lid for superior inventory control. Two and three inch high boxes are available and accommodate optional fiberboard dividers listed below.



82007-142, 82007-144, 82007-162

Fiberboard Dividers (1" high to fit 5" x 5" boxes)

VWR Cat#	Style	Dimensions (In.)	List Pricing
82007-156	100 cell	5" x 5" x 1"	\$1.40/ ea.
82007-154	81 cell	5" x 5" x 1"	\$1.40/ ea.
82007-152	64 cell	5" x 5" x 1"	\$1.40/ ea.
82007-150	49 cell	5" x 5" x 1"	\$1.40/ ea.
82007-148	25 cell	5" x 5" x 1"	\$1.40/ ea.
82007-146	16 cell	5" x 5" x 1"	\$1.40/ ea.

Cryogenic fiberboard dividers can be utilized in any 5"x5" mechanical or liquid nitrogen cryogenic freezer box. Fiberboard dividers are available in configurations from 16 to 100 cells.

VWR Family Brands: Comprehensive product offerings designed to meet a wide variety of research needs.



VWR™
INTERNATIONAL

vwr.com 1.800.932.5000

Order from Over 750,000 Products

Sales & Inventory Locations

Pacific Northwest

Anchorage, AK
Boise, ID
Brisbane, CA
Seattle, WA
Tualatin, OR

Southwest

Albuquerque, NM
Aurora, CO
Phoenix, AZ
San Diego, CA
San Dimas, CA

Midwest

Batavia, IL
Indianapolis, IN
Minneapolis, MN
St. Louis, MO

Gulf

Lake Charles, LA
Sugarland, TX

Northeast

Franklin, MA
Rochester, NY

Mid-Atlantic

Bridgeport, NJ
Clarksburg, MD

Southeast

Morrisville, NC
Suwanee, GA

Canada

Mississauga, Ontario
Ville Mont-Royal, Québec

Call Direct for Specialized Service Locations

International

911 Commerce Ct.
Buffalo Grove, IL 60089
847-520-6170
fax: 847-229-0788

Puerto Rico

Carr. #869, Km. 1.5 M-4
Royal Industrial Park
Catano, PR 00962
800-932-5000
787-275-4500
fax: 787-788-4320

Global and

Corporate Accounts

1310 Goshen Pkwy.
West Chester, PA 19380
610-431-1700
fax: 610-431-9174

Mexico

Super Av. Lomas Verdes
No. 464 2o Piso Col.
Los Alamos C.P. 53230
Naucalpan de Juarez,
Edo. de Mexico 2122-1600
011-52-552-122-1630
fax: 011-52-555-344-7532

Production Supplies and Services

975 Overland Ct.
San Dimas, CA 91773
877-VWR-PSCS
(877-897-7727)
fax: 888-890-9124

Furniture Group

3021 Gateway Drive
Suite 280
Irving, TX 75063
888-624-2432
fax: 972-753-1389

VWR Laboratory Distributors Services

911 Commerce Ct.
Buffalo Grove, IL 60089
800-873-8971
fax: 847-229-0668

Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.

©2004 VWR International, Inc. All rights reserved. Printed in U.S.A.

0304(10M) Lit. No. 99991