



VWR™ Disposable Pellet Mixers

Lightweight
Battery Powered
Disposable



Pestle Motor Mixer - 47747-370
1.5 ml Micro-tube & Pestle - 47747-366

VWR Disposable Pellet Mixer

- Lightweight
- Battery Powered
- Disposable

Designed to homogenize cells and plant tissue or to re-suspend protein and DNA pellets. Manufactured from PBTB in class 10,000 clean room, VWR pestles and tubes are DNAase – RNAase and Pyrogen free and individually wrapped to guarantee purity. Their grooved grip handle is tailored to fit with our micro tube. It is easy to turn and does not slip even when wearing gloves. Autoclavable: Pestle Length: 7.5cm

The VWR Pellet Mixer is a lightweight battery powered instrument that can be used with disposable pestles, mixers and micro tubes. Cell tissue and pellets are completely homogenized within seconds by high speed vortexing action that greatly reduces hand fatigue. Comes complete with 2 "AA" batteries and Pestle/Mixer adapter.

Disposable Pellet Mixer Specifications

Catalog#	Description		
47747-358	1.5 ml Pestle	PK/100	\$71.95
47747-360	.5 ml Pestle	PK/100	\$71.95
47747-362	1.5 ml Micro-tube	PK/100	\$21.50
47747-364	.5ml Micro-tube	PK/100	\$21.50
47747-366	1.5 ml Micro-tube & Pestle	PK/100	\$88.50
47747-368	.5ml Micro-tube & Pestle	PK/100	\$88.50
47747-370	Pestle Motor Mixer	1/EA	\$79.36
47747-372	Motor Adapter Replacement	1/EA	\$6.00



VWR  TM
INTERNATIONAL

vwr.com 1.800.932.5000

Order from Over 750,000 Products

Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.

©2003 VWR International, Inc. All rights reserved. Printed in U.S.A.



VWR Products:
Providing value at competitive prices