



## VWR Galaxy Digital Microcentrifuges

**The Galaxy**

**has**

**Gone**

**Digital**



*37001-594, 37001-300*

**Galaxy**

## VWR Galaxy Digital Microcentrifuge

- Digital control panel
- Compact design
- Five speed ranges

The newest additions to the VWR Galaxy Microcentrifuge line are five digital models. A variety of speed ranges, from 5,000 to 16,000 x g are available to meet all your microcentrifuge needs.

Large digital displays and membrane keypads allow for precise setting of operating parameters and accurate representation of their values during a run. Speed/rcf can be set and displayed in rpm or g.

In addition to digital control, the larger Galaxy models (14D, 16D and 16DH) feature brushless, maintenance-free motors for faster acceleration and quieter performance. The entire line has been designed to require a minimum of benchspace. Each of the VWR Galaxy Digital Microcentrifuges is supplied complete with a 1.5/2.0ml rotor.

### Specifications

	5D	7D	14D	16D	16DH
Max. speed	8,500rpm	10,000rpm	13,000rpm	14,000rpm	13,300rpm
Max. RCF	5,200 x g	7,200 x g	14,000 x g	16,000 x g	16,000 x g
Rotor	12 place	12 place	18 place	18 place	24 place
Timer	15 min	99 min	99 min	99 min	99 min
Dimensions	7.5 x 7 x 4.5" 191 x 178 x 127 mm		10.5 x 9.3 x 7" 267 x 236 x 178 mm		
Weight	1.5 kg	1.5 kg	5.5 kg	5.5 kg	5.9 kg
Electrical	120V, 50-60 Hz				

### Ordering Information

Cat. No.	Description	Price
37001-594	Galaxy 5D with 12 x 1.5/2.0ml rotor	\$695.00
37001-294	Galaxy 7D with 12 x 1.5/2.0ml rotor	895.00
37001-296	Galaxy 14D with 18 x 1.5/2.0ml rotor	1350.00
37001-298	Galaxy 16D with 18 x 1.5/2.0ml rotor	1795.00
37001-300	Galaxy 16DH with 24 x 1.5/2.0ml rotor	1895.00
20668-224	Adapters for 0.5ml tubes, pack of 6	35.00
20668-226	Adapters for 0.4ml tubes, pack of 6	35.00
20668-228	Adapters for 0.2ml tubes, pack of 6	35.00

The economical Galaxy 7 and Galaxy 16 are also available from VWR International

VWR Family Brands: Comprehensive product offerings designed to meet a wide variety of research needs.



**vwr.com** 1.800.932.5000

Order from Over 400,000 Products

### Sales & Inventory Locations

#### Pacific Northwest

Anchorage, AK  
Boise, ID  
Portland, OR  
San Francisco, CA  
Seattle, WA

#### Southwest

Albuquerque, NM  
Denver, CO  
Phoenix, AZ  
San Diego, CA  
San Dimas, CA

#### Midwest

Chicago, IL  
Indianapolis, IN  
Minneapolis, MN  
St. Louis, MO

#### Gulf States

Dallas, TX  
Houston, TX  
Lake Charles, LA

#### Northeast

Boston, MA  
Cleveland, OH  
New York, NY  
Pittsburgh, PA  
Rochester, NY

#### Mid-Atlantic

Baltimore, MD  
Bridgeport, NJ  
Columbus, OH

#### Southeast

Atlanta, GA  
Miami, FL  
Research Triangle Park, NC  
Oak Ridge, TN

#### VWR Canada

Mississauga, Ontario  
Ville Mont-Royal, Québec

### Call Direct for Specialized Service Locations

#### International

911 Commerce Ct.  
Buffalo Grove, IL 60089  
847-520-6170  
fax: 847-229-0788

#### Puerto Rico

Carr. #869, Km. 1.5 M-4  
Royal Industrial Park  
Catano, PR 00962  
800-932-5000  
787-275-4500  
fax: 787-788-4320

#### Switzerland

Ruchligstrasse 20  
P.O. Box 964, CH-8953  
Dietikon, Switzerland  
011-41-1-745-1155  
fax: 011-41-1-745-1150

#### Corporate Accounts

1310 Goshen Pkwy.  
West Chester, PA 19380  
610-431-1700  
fax: 610-431-9174

#### Production Supplies and Services

975 Overland Ct.  
San Dimas, CA 91773  
877-VWR-PSCS  
(877-897-7727)  
fax: 888-890-9124

#### Furniture Division

P.O. Box 3405  
Irving, TX 75015  
888-624-2432  
fax: 972-753-1389

#### VWR Direct

911 Commerce Ct.  
Buffalo Grove, IL 60089  
800-444-0880  
fax: 847-229-8465

#### VWR Laboratory Distributors Services

911 Commerce Ct.  
Buffalo Grove, IL 60089  
800-873-8971  
fax: 847-229-0668

Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.