



VWR™ Heavy Duty Support Jacks

Ultra stable support
when and where
it's needed



(Clockwise) 14233-372, 14233-368 w/14217-532, 14233-232 14233-360, 14233-364

VWR™ Heavy Duty Support Jacks

- Exceptional stability and durability
- Stainless steel construction
- Seven convenient sizes to choose from
- Autoclavable and chemical resistant
- Optional support rod and ratcheting handle

These heavy duty, stainless steel support jacks are ultra stable lifting platforms with exceptional strength and durability. Constructed of 300 series stainless steel, the support jacks are designed for use in extreme environments and high load applications. The support jacks' durable construction allows them to be autoclaved or chemically cleaned. Support jacks are equipped with over-sized, sure-grip adjustment knobs that provide extra leverage and easy height adjustments.

These lifts are ideal for use in fume hoods or bench tops and hold a variety of items such as glassware, hotplates, baths, and magnetic stirrers. VWR offers seven different sizes to accommodate any application.

Optional Accessories

Support Rod Kit- This kit allows you to mount a 17" vertical support rod to the upper deck of a 6x6", 8x8", 10x10", 12x12", and 16x16" support jacks. Ideal for creating an adjustable support stand or mounting various items such as thermometer clamps, temperature probes, flask and column clamps.

Ratcheting handle tool- Designed to add extra leverage to your 12x12" and 16x16" support jacks. This recommended ratcheting handle, easily attaches to the actuating rod to allow easy, accurate adjustments.

Specifications

| Deck Size in(cm) | Max Ht. in(cm) | Closed Ht in(cm) | Max Load* lb(kg) |
|-------------------------|-------------------|---------------------|---------------------|
| 3 x 3 (7.62 x 7.62) | 5.00 (12.70) | 2.50 (6.35) | 20 (9.07) |
| 4 x 4 (10.16 x 10.16) | 5.00 (12.70) | 2.50 (6.35) | 30 (13.61) |
| 6 x 6 (15.24 x 15.24) | 9.75 (24.76) | 3.00 (7.62) | 35(15.88) |
| 8 x 8 (20.32 x 20.32) | 13.00 (33.02) | 3.50 (8.89) | 55 (24.95) |
| 10 x 10 (25.40 x 25.40) | 13.00 (33.02) | 3.50 (8.89) | 80 (36.29) |
| 12 x 12 (30.48 x 30.48) | 19.50 (49.53) | 4.00 (10.16) | 150 (68.04) |
| 16 x 16 (40.64 x 40.64) | 19.50 (49.53) | 4.00 (10.16) | 180 (81.65) |

17" Support rod kit - includes 17" rod, 1-jam nut, 1-flat washer

Ratchet handle tool - for use on 12x12 and 16x16 support jacks

* Max load rating represents static weight only. Static weight is the amount a unit can hold, not lift.

| Cat. No. | Description | Price |
|-----------|----------------------|-----------|
| 14233-360 | 3 x 3 support jack | \$ 136.00 |
| 14233-362 | 4 x 4 support jack | 146.00 |
| 14233-364 | 6 x 6 support jack | 219.00 |
| 14233-366 | 8 x 8 support jack | 278.00 |
| 14233-368 | 10 x 10 support jack | 315.00 |
| 14233-370 | 12 x 12 support jack | 607.00 |
| 14233-372 | 16 x 16 support jack | 1099.00 |
| 14217-532 | 17" Support rod kit | 19.99 |
| 14233-232 | Ratchet tool | 39.99 |

VWR:

Providing value at competitive prices

VWR
INTERNATIONAL

vwr.com 1.800.932.5000

Order from Over 750,000 Products

Sales & Inventory Locations

Pacific Northwest

Anchorage, AK
Boise, ID
Brisbane, CA
Seattle, WA
Tualatin, OR

Southwest

Albuquerque, NM
Aurora, CO
Phoenix, AZ
San Diego, CA
San Dimas, CA

Midwest

Batavia, IL
Indianapolis, IN
Minneapolis, MN
St. Louis, MO

Gulf

Houston, TX
Lake Charles, LA
Sugarland, TX

Northeast

Marlborough, MA
Rochester, NY

Mid-Atlantic

Bridgeport, NJ
Clarksburg, MD

Southeast

Miami, FL
Morrisville, NC
Oak Ridge, TN
Suwanee, GA

Canada

Mississauga, Ontario
Ville Mont-Royal,
Québec

Call Direct for Specialized Service Locations

International

911 Commerce Ct.
Buffalo Grove, IL 60089
847-520-6170
fax: 847-229-0788

Puerto Rico

Carr. #869, Km. 1.5 M-4
Royal Industrial Park
Catano, PR 00962
800-932-5000
787-275-4500
fax: 787-788-4320

Global and Corporate Accounts

1310 Goshen Pkwy.
West Chester, PA 19380
610-431-1700
fax: 610-431-9174

Production Supplies and Services

975 Overland Ct.
San Dimas, CA 91773
877-VWR-PSCS
(877-897-7727)
fax: 888-890-9124

Furniture Division

3021 Gateway Drive
Suite 280
Irving, TX 75063
888-624-2432
fax: 972-753-1389

VWR Laboratory Distributors Services

911 Commerce Ct.
Buffalo Grove, IL 60089
800-873-8971
fax: 847-229-0668

Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.

©2003 VWR International, Inc. All rights reserved. Printed in U.S.A.

01/03

