



VWR™ Lab Disc

Now

Thinner


Than

Ever!!!



#82007-828

So thin it fits in your pocket!



Actual Size - #82007-828

VWR™ Lab Disc

- No moving parts
- Rotation direction is reversible
- Weighs less than 1 lb.
- Explosion hazard zone 2

VWR™ Lab Disc, the innovative, ultra-flat magnetic stirrer from VWR™ International. Advanced magnetic coil technology with no moving parts makes the Lab Disc drive completely wear-free. If desired, the Lab Disc can reverse rotation direction automatically every 30 seconds to ensure better mixing. A reliable, strong mixer with a great design to bring color into your lab. The Lab Disc is perfect for all of your small stirring applications.

Stirring quantity: 800 mL
Power input: 3W
Speed range: 15 - 1,500 rpm
Material: polyester
Dimensions: 114 x 160 x 10 mm
Weight: 0.3 kg

Cat. No.
82007-828

List Price
\$ 154.00

VWR Products:
Providing value at competitive prices

VWR 
INTERNATIONAL

vwr.com 1.800.932.5000

Order from Over 750,000 Products

Sales & Inventory Locations

Pacific Northwest

Anchorage, AK
Boise, ID
Brisbane, CA
Seattle, WA
Tualatin, OR

Southwest

Albuquerque, NM
Aurora, CO
Phoenix, AZ
San Diego, CA
San Dimas, CA

Midwest

Batavia, IL
Indianapolis, IN
Minneapolis, MN
St. Louis, MO

Gulf

Houston, TX
Lake Charles, LA
Sugarland, TX

Northeast

Marlborough, MA
Rochester, NY

Mid-Atlantic

Bridgeport, NJ
Clarksburg, MD

Southeast

Morrisville, NC
Suwanee, GA

Canada

Mississauga, Ontario
Ville Mont-Royal,
Québec

Call Direct for Specialized Service Locations

International

911 Commerce Ct.
Buffalo Grove, IL 60089
847-520-6170
fax: 847-229-0788

Puerto Rico

Carr. #869, Km. 1.5 M-4
Royal Industrial Park
Catano, PR 00962
800-932-5000
787-275-4500
fax: 787-788-4320

Global and Corporate Accounts

1310 Goshen Pkwy.
West Chester, PA 19380
610-431-1700
fax: 610-431-9174

Mexico

Super Av. Lomas Verdes
No. 464 2o Piso Col.
Los Alamos C.P. 53230
Naucalpan de Juarez,
Edo, de Mexico 2122-1600
011-52-552-122-1630
fax: 011-52-555-344-7532

Production Supplies and Services

975 Overland Ct.
San Dimas, CA 91773
877-VWR-PSCS
(877-897-7727)
fax: 888-890-9124

Furniture Group

3021 Gateway Drive
Suite 280
Irving, TX 75063
888-624-2432
fax: 972-753-1389

VWR Laboratory Distributors Services

911 Commerce Ct.
Buffalo Grove, IL 60089
800-873-8971
fax: 847-229-0668

Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.

©2003 VWR International, Inc. All rights reserved. Printed in U.S.A.

10/03

11/03 (10m) Lit. No. 99967

