



VWR™ Nutating Mixer

A new
dimension in
mixing and
shaking



82007-202

VWR™ Nutating Mixer

- Three dimensional motion
- Provides thorough yet gentle mixing
- Ideal for blood tubes and centrifuge tubes

The VWR Nutating Mixer combines the actions of an orbital shaker and a rocker to produce a three dimensional motion. This resulting motion provides thorough mixing, yet is gentle enough for fragile samples.

The speed and tilt angle of the nutating mixer are fixed to provide the ideal conditions for mixing and resuspending samples in blood tubes and centrifuge tubes.

A dimpled mat is supplied with the VWR Nutating Mixer to hold tubes in place. The mat will hold as many as 21 blood collection tubes, 30 pediatric blood collection tubes, 60 microtubes, seventeen 15mL centrifuge tubes or nine 50mL centrifuge tubes.

The VWR Nutating Mixer has a small footprint and low profile, requiring a minimum amount of bench space. For temperature sensitive samples, the mixer may be placed in a cold room.

Specifications

Speed range	24rpm*
Motion	Three dimensional
Tilt range	Fixed, 20°
Maximum load	0.8 kg/1.75 lb
Ambient operating range	+4° to 40°C
Platform size	24 x 17.15 cm/9.5 x 6.75 inches
Dimensions (WxDxH)	24 x 17.15 x 15/24 cm/9.75 x 6.75 x 6 inches
Weight	.88 kg/1.9 lb
Electrical	120V, 60Hz

*Max. speed for 230V version is 20rpm.

Ordering Information

Cat. No.	Description	Price
82007-202	VWR Nutating Mixer with dimpled mat, 120V, 60 Hz	\$425.00
83007-210	VWR Nutating Mixer with dimpled mat, 230V, 50 Hz	425.00



vwr.com 1.800.932.5000

Order from Over 750,000 Products

Sales & Inventory Locations

Pacific Northwest Anchorage, AK Boise, ID Brisbane, CA Seattle, WA Tualatin, OR	Gulf Houston, TX Lake Charles, LA Sugarland, TX
Southwest Albuquerque, NM Aurora, CO Phoenix, AZ San Diego, CA San Dimas, CA	Northeast Marlborough, MA Rochester, NY
Midwest Batavia, IL Indianapolis, IN Minneapolis, MN St. Louis, MO	Mid-Atlantic Bridgeport, NJ Clarksburg, MD
	Southeast Morrisville, NC Suwanee, GA
	Canada Mississauga, Ontario Ville Mont-Royal, Québec

Call Direct for Specialized Service Locations

International 911 Commerce Ct. Buffalo Grove, IL 60089 847-520-6170 fax: 847-229-0788	Production Supplies and Services 975 Overland Ct. San Dimas, CA 91773 877-VWR-PSCS (877-897-7727) fax: 888-890-9124
Puerto Rico Carr. #869, Km. 1.5 M-4 Royal Industrial Park Catano, PR 00962 800-932-5000 787-275-4500 fax: 787-788-4320	Furniture Group 3021 Gateway Drive Suite 280 Irving, TX 75063 888-624-2432 fax: 972-753-1389
Global and Corporate Accounts 1310 Goshen Pkwy. West Chester, PA 19380 610-431-1700 fax: 610-431-9174	VWR Laboratory Distributors Services 911 Commerce Ct. Buffalo Grove, IL 60089 800-873-8971 fax: 847-229-0668
Mexico Super Av. Lomas Verdes No. 464 2o Piso Col. Los Alamos C.P. 53230 Naucaipan de Juarez, Edo, de Mexico 2122-1600 011-52-552-122-1630 fax: 011-52-555-344-7532	

Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.

©2003 VWR International, Inc. All rights reserved. Printed in U.S.A.

10/03

11/03 (5M) Lit. No. 99945

VWR:
Providing value at competitive prices.

