

More
Than Meets
the Eye



12777-984

VWR VistaVision™ Refractometers

The VistaVision™ line of refractometers provides accurate readings of juices, beverages, industrial fluids, and salt water, and is designed for veterinary and automotive applications. Available with ATC (Automatic Temperature Compensation), which is ideal for field use.

The Brix Series offers a direct reading of the Brix scale. This scale is used for sugar-related liquids, such as fruit, soft drinks, and wine. This unit can also be used for finding out the concentration of industrial fluids, such as automotive coolant.

VWR Catalog Number	Description	Divisions	Accuracy ±	List Price
12777-978	Refractometer Brix 0-10 - ATC	0.1°	0.1°	\$225.00
12777-980	Refractometer Brix 0-18	0.2°	0.15°	\$175.00
12777-982	Refractometer Brix 0-32	0.2°	0.2°	\$175.00
12777-984	Refractometer Brix 0-32 - ATC	0.2°	0.2°	\$225.00
12777-986	Refractometer Brix 28-62	0.2°	0.2°	\$195.00
12777-988	Refractometer Brix 45-82	0.5°	0.3°	\$205.00
12777-990	Refractometer Brix / Honey 58-90	1.0°	0.5°	\$205.00

The Automotive Series is designed for testing antifreeze or battery fluids. This model has a viewing area with three scales for the different fluids: Propylene Glycol, Ethylene Glycol and Battery Charge level. This model tests the freezing point of antifreeze or the recharge level of battery fluid.

VWR Catalog Number	Description	Divisions	Accuracy ±	List Price
12778-000	Refractometer Automotive			\$225.00
	- Propylene Glycol	10°	±5.0°	
	- Ethylene Glycol	10°	±5.0°	
	- Battery Recharge Level	n/a	n/a	

The Salt Series provides a direct reading of the specific gravity and concentration (ppt) of salt water. This unit is designed for food processing applications.

VWR Catalog Number	Description	Divisions	Accuracy ±	List Price
12777-992	Refractometer Salt --100ppt - ATC			\$205.00
	- 0-100ppt Salinity	1.0 ppt	±1.0 ppt	
	- 1.0 to 1.07 Specific Gravity	0.2%	0.2%	
12777-994	Refractometer Salt 0-28	0.2%	0.2%	\$180.00

The Veterinary Series is intended for veterinary use only. These units offer veterinarians a quick indication of vital fluid levels. This unit offers three readings of urine, total serum, and refractive index.

VWR Catalog Number	Description	Divisions	Accuracy ±	List Price
12777-996	Refractometer Veterinary			\$170.00
	- Urine Specific Gravity	0.002	±0.002	
	- Serum Protein	0.2 g/dl	±0.2 g/dl	
	- Refractive Index	0.0005	±0.002	
12777-998	Refractometer Veterinary - ATC			\$170.00
	- Urine Specific Gravity	0.002	±0.002	
	- Serum Protein	0.2 g/dl	±0.2 g/dl	
	- Refractive Index	0.0005	±0.002	



vwr.com 1.800.932.5000

Order from Over 400,000 Products

Sales & Inventory Locations

Pacific Northwest

Anchorage, AK
Boise, ID
Brisbane, CA
Seattle, WA
Tualatin, OR

Southwest

Albuquerque, NM
Aurora, CO
Phoenix, AZ
San Diego, CA
San Dimas, CA

Midwest

Batavia, IL
Indianapolis, IN
Minneapolis, MN
St. Louis, MO

Gulf

Houston, TX
Lake Charles, LA
Sugarland, TX

Northeast

Marlborough, MA
Rochester, NY

Mid-Atlantic

Bridgeport, NJ
Clarksburg, MD

Southeast

Miami, FL
Morrisville, NC
Oak Ridge, TN
Suwanee, GA

Canada

Mississauga, Ontario
Ville Mont-Royal,
Québec

Call Direct for Specialized Service Locations

International

911 Commerce Ct.
Buffalo Grove, IL 60089
847-520-6170
fax: 847-229-0788

Puerto Rico

Carr. #869, Km. 1.5 M-4
Royal Industrial Park
Catano, PR 00962
800-932-5000
787-275-4500
fax: 787-788-4320

Global and Corporate Accounts

1310 Goshen Pkwy.
West Chester, PA 19380
610-431-1700
fax: 610-431-9174

Production Supplies and Services

975 Overland Ct.
San Dimas, CA 91773
877-VWR-PSCS
(877-897-7727)
fax: 888-890-9124

Furniture Division

3021 Gateway Drive
Suite 280
Irving, TX 75063
888-624-2432
fax: 972-753-1389

VWR Laboratory Distributors Services

911 Commerce Ct.
Buffalo Grove, IL 60089
800-873-8971
fax: 847-229-0668

Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.

©2002 VWR International, Inc. All rights reserved. Printed in U.S.A.

05/02

0602 (7.5M) Lit. No. 99858

VWR Family Brands: Comprehensive product offerings designed to meet a wide variety of research needs.

