



VWR™ Rocking Platform

Maximum
efficiency in
staining and
blot washing



40000-304

VWR™ Rocking Platforms

- Gentle rocking effect lifts and washes over samples
- Durable construction with stainless steel platforms
- Double platform model for increased capacity

The VWR Rocking Platforms feature variable speed and adjustable tilt to provide optimal control over rocking conditions. Parameters can be set to provide the gentle motion required for staining fragile gels as well as the more vigorous conditions required for washing membranes. Should an application require temperatures above ambient, the rocking platforms can be used in temperature controlled environments, up to 65°C.

A heavy duty motor and stainless steel rocking platforms can support loads of up to 7 kilograms. The closed-loop, solid state speed control ensures the rocker operates at constant speed, regardless of voltage or load fluctuations. Permanently lubricated ball bearings provide years of maintenance-free operation.

For added capacity, the VWR Rocking Platform is available with two stacked platforms. Clearance between the platforms is 9 cm. Both the single and double platform models are supplied with nonslip rubber mats to keep sample containers in place.

Specifications

Speed range	4 to 160 tilts/minute*
Motion	Rocking
Tilt range	Adjustable, ±20 °
Clearance between platforms	9 cm/3.5 inches
Maximum load	7 kg/15.4 lb
Ambient operating range	+4° to 65°C
Dimensions (WxDxH)	35.6 x 33 x 16.5 cm/14 x 13 x 6.5 inches
Weight	5 kg/11 lb
Electrical	120V or 230V, 50-60Hz

*Maximum speed for 230V unit is 150 tilts/minute.

Ordering Information

Cat. No.	Description	Price
40000-300	VWR Model 100 single rocking platform, 120V, 50/60 Hz	\$720.00
40000-302	VWR Model 100 single rocking platform, 230V, 50/60 Hz	\$775.00
40000-304	VWR Model 200 double rocking platform, 120V, 50/60 Hz	\$800.00
40000-306	VWR Model 200 double rocking platform, 230V, 50/60 Hz	\$850.00



vwr.com 1.800.932.5000

Order from Over 750,000 Products

Sales & Inventory Locations

Pacific Northwest Anchorage, AK Boise, ID Brisbane, CA Seattle, WA Tualatin, OR	Gulf Houston, TX Lake Charles, LA Sugarland, TX
Southwest Albuquerque, NM Aurora, CO Phoenix, AZ San Diego, CA San Dimas, CA	Northeast Marlborough, MA Rochester, NY
Midwest Batavia, IL Indianapolis, IN Minneapolis, MN St. Louis, MO	Mid-Atlantic Bridgeport, NJ Clarksburg, MD
	Southeast Miami, FL Morrisville, NC Oak Ridge, TN Suwanee, GA
	Canada Mississauga, Ontario Ville Mont-Royal, Québec

Call Direct for Specialized Service Locations

International 911 Commerce Ct. Buffalo Grove, IL 60089 847-520-6170 fax: 847-229-0788	Production Supplies and Services 975 Overland Ct. San Dimas, CA 91773 877-VWR-PSCS (877-897-7727) fax: 888-890-9124
Puerto Rico Carr. #869, Km. 1.5 M-4 Royal Industrial Park Catano, PR 00962 800-932-5000 787-275-4500 fax: 787-788-4320	Furniture Division 3021 Gateway Drive Suite 280 Irving, TX 75063 888-624-2432 fax: 972-753-1389
Global and Corporate Accounts 1310 Goshen Pkwy. West Chester, PA 19380 610-431-1700 fax: 610-431-9174	VWR Laboratory Distributors Services 911 Commerce Ct. Buffalo Grove, IL 60089 800-873-8971 fax: 847-229-0668

VWR:

Providing value at competitive prices.



Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.

©2003 VWR International, Inc. All rights reserved. Printed in U.S.A.

01/03

06/03 (5M) Lit. No. 99954