



VWR™ Safe-T-Flex™ Caps

The Ultimate

Closure for

Personal Safety and

Sample Integrity



VWR Cat. Nos.: 60828-762 and 60828-768

VWR™ Safe-T-Flex™ Caps

- Will not aerosol when opened.
- Fits all varieties of blood collecting and disposable culture tubes.
- Liquid-tight seal prevents sample evaporation or spillage.
- Withstands centrifugation, agitation, refrigeration and freezing.
- Easy, one-hand removal and application.
- Available in seven colors. White colored caps are trace-metal free.

Now there is a safe, economical way to cap and recap blood collection tubes and disposable glass culture tubes, regardless of manufacturer. VWR Safe-T-Flex caps feature a unique, flexible, over-locking design that is far superior to Parafilm*. Available for a wide range of biomedical and industrial applications, VWR Safe-T-Flex caps lock securely on the outside of the tube lip and can be removed with just one hand.

The VWR Safe-T-Flex liquid-tight seal ensures sample integrity and prevents spillage of serum samples. This special seal also helps to protect your laboratory staff from the aerosols of highly infectious organisms such as TB and the HTLV-III virus. VWR Safe-T-Flex caps are available in seven colors, which helps guard against sample cross-contamination, and is convenient for color-coding your particular research and development projects. You can also write directly on the VWR Safe-T-Flex cap with indelible ink.

VWR Safe-T-Flex caps are packaged in bags of 1,000 and a convenient dispenser is also available for storage and easy access.

For 10mm Blood Collection and Culture Tubes

VWR Cat. No.	Color	List Price/1,000
60828-700	White	\$49.16
60828-702	Red	49.16
60828-704	Lavender	49.16
60828-706	Green	49.16
60828-708	Blue	49.16
60828-710	Grey	49.16
60828-712	Yellow	49.16

For 13mm Culture Tubes

VWR Cat. No.	Color	List Price/1,000
60828-738	White	\$49.16
60828-740	Red	49.16
60828-742	Lavender	49.16
60828-744	Green	49.16
60828-746	Blue	49.16
60828-748	Grey	49.16
60828-750	Yellow	49.16

For 12mm Culture Tubes/13mm Blood Collection Tubes

VWR Cat. No.	Color	List Price/1,000	10pk./cs.
60828-720	White	\$49.16	—
60828-722	Red	49.16	—
60828-724	Lavender	49.16	\$491.60
60828-726	Green	49.16	\$491.60
60828-728	Blue	49.16	\$491.60
60828-730	Grey	49.16	—
60828-732	Yellow	49.16	—

For 16mm Blood Collection & Culture Tubes

VWR Cat. No.	Color	List Price/1,000	10pk./cs.
60828-760	White	\$49.16	—
60828-762	Red	49.16	\$491.60
60828-764	Lavender	49.16	\$491.60
60828-766	Green	49.16	—
60828-768	Blue	49.16	\$491.60
60828-770	Grey	49.16	—
60828-772	Yellow	49.16	—



Using a VWR Safe-T-Flex cap is an easy, one-hand operation. To close, press down firmly with thumb to lock in place.

To open, simply pry cap up with thumb.



*Parafilm is a registered trademark of American Can. Co./Dixie Marathon.

VWR:
Providing value at competitive prices



Order from Over 400,000 Products

Sales & Inventory Locations

Pacific Northwest

Anchorage, AK
Boise, ID
Portland, OR
San Francisco, CA
Seattle, WA

Southwest

Albuquerque, NM
Denver, CO
Phoenix, AZ
San Diego, CA
San Dimas, CA

Midwest

Chicago, IL
Indianapolis, IN
Minneapolis, MN
St. Louis, MO

Gulf States

Dallas, TX
Houston, TX
Lake Charles, LA

Northeast

Boston, MA
Cleveland, OH
Marlboro, MA
New York, NY
Pittsburgh, PA
Rochester, NY

Mid-Atlantic

Baltimore, MD
Bridgeport, NJ
Columbus, OH

Southeast

Atlanta, GA
Miami, FL
Research Triangle Park, NC
Oak Ridge, TN

VWR Canada

Mississauga, Ontario
Ville Mont-Royal, Québec

Call Direct for Specialized Service Locations

International

911 Commerce Ct.
Buffalo Grove, IL 60089
847-520-6170
fax: 847-229-0788

Puerto Rico

Carr. #869, Km. 1.5 M-4
Royal Industrial Park
Catano, PR 00962
800-932-5000
787-275-4500
fax: 787-788-4320

Switzerland

Ruchligstrasse 20
P.O. Box 964, CH-8953
Dietikon, Switzerland
011-41-1-745-1155
fax: 011-41-1-745-1150

Corporate Accounts

1310 Goshen Pkwy.
West Chester, PA 19380
610-431-1700
fax: 610-431-9174

Production Supplies and Services

975 Overland Ct.
San Dimas, CA 91773
877-VWR-PSCS
(877-897-7727)
fax: 888-890-9124

Furniture Division

P.O. Box 3405
Irving, TX 75015
888-624-2432
fax: 972-753-1389

VWR Direct

911 Commerce Ct.
Buffalo Grove, IL 60089
800-444-0880
fax: 847-229-8465

VWR Laboratory Distributors Services

911 Commerce Ct.
Buffalo Grove, IL 60089
800-873-8971
fax: 847-229-0668

Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.

Safe-T-Flex Scientific Holdings Corp.
©2002 VWR International, Inc. All rights reserved. Printed in U.S.A.

1/02

0402 (10M) Lit. No. 99861