

VWR™ Traceable® IR Thermometers

Take temperatures without ever touching your sample— Easy to use IR thermometers, just point and read the temperature. Samples remain contamination-free and sterile solutions remain sterile with probeless measurements. It will even take temperatures of moving objects.

VWR Cat #	21800-018	77776-724	12777-846	12777-844	33501-416
Traceable® NIST Cert Supplied	yes	yes	yes	yes	yes
Range	14 to 572°F -10 to 300°C	14 to 662°F -10 to 350°C	-4 to 788°F -20 to 420°C	-0 to 600°F -20 to 315°C	-7 to 230°F -22 to 110°C
Resolution	0.1/ 1°	0.1°	1°	1°	0.1°
Accuracy (0 to 100°C)	±2.0°C	±2.0°C	±2°C	±2°C	±1.0°C
Unit Size (inch)	7x3x1.25	7.75x4.75x2.3	7x2x1.5	6x2x1.25	2.5x1.5x0.5
Alarm	no	no	no	yes	no
Min/Max	yes	yes	no	yes	no
List Price	\$269.66	\$146.95	\$103.95	\$116.73	\$34.95



VWR Cat. No. **33501-416**
Traceable® Micro IR Thermometer Pocket size, clock feature, simple operation, take anywhere—



VWR Cat. No. **21800-018**
Traceable® Digital Thermometer Probe with cable allows placing sensor anywhere, RS232 serial output, ultra-fast sampling—



VWR Cat. No. **77776-724**
Traceable® Noncontact Thermometer Features adjustable emissivity, data hold, min/max memory and accessory Type-K probe—



VWR Cat. No. **12777-846**
Traceable® Infrared Gun Thermometer Laser-beam sightings permits one-handed “point & shoot” readings, hold and backlight—



VWR Cat. No. **12777-844**
Traceable® Memory/Alarm Thermometer Backlight, 9-point memory, laser targeting, high/low alarm—

All thermometers are supplied with a Traceable® Certificate, traceable to NIST (National Institute of Standards and Technology). It meets the exacting requirements of ISO/IEC 17025. Calibration supplied from an ISO 9001 Quality Certified Manufacturer and an ISO 17025 A2LA Accredited Calibration Laboratory.

VWR:
Providing value at competitive prices

VWR™
INTERNATIONAL

vwr.com 1.800.932.5000

Order from Over 750,000 Products

Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.


©2003 VWR International, Inc. All rights reserved. Printed in U.S.A.





VWR™ Traceable® IR Thermometers

The newest
way to measure
temperatures--
infrared



VWR Cat. No. 33501-416