



## VWR™ pH Calibration Buffers

For pH Measurements  
That Are Important  
To Your Work



## VWR pH Calibration Buffers

- **Wide Product Range** - pH 1.00 - pH 12.00
- **Quality Assured** - Accuracy + 0.02 pH, NIST traceable, Expiration Dated
- **Documentation** - Certificates of Analysis Material Safety Data Sheets

VWR No.	Description	Size	List Price Each	Case (6 per cs.)
34170-192	Multipak 2 each of pH 4,7,10	6 x 500ml	\$92.74	
34170-199	Multipak pH 4,7,10	20 x 20ml	37.54	
34170-218	Buffer pH 1.00	500ml	8.11	
34170-222		20L	114.41	
34170-182	Buffer pH 1.68	500ml	29.39	\$176.34
34170-100	Buffer pH 2.00	500ml	27.76	165.84
34170-228		20L	126.01	
34170-103	Buffer pH 3.00	500ml	27.64	165.84
34170-230		20L	147.00	
34170-106	Buffer pH 4.00	500ml	19.65	117.90
34170-164		20L	87.55	
34170-127	Buffer pH 4.00 Red	500ml	19.65	117.90
34170-155		20L	87.50	
34170-172		20 x 20ml	60.47	
34170-109	Buffer pH 5.00	500ml	27.64	165.84
34170-210		20l	99.92	
34170-112	Buffer pH 6.00	500ml	26.55	159.30
34170-115	Buffer pH 7.00	500ml	18.67	112.02
34170-167		20L	87.46	
34170-130	Buffer pH 7.00 Yellow	500ml	19.65	117.90
34170-158		20L	87.47	
34170-175		20 x 20ml	58.77	
34170-118	Buffer pH 8.00	500ml	26.55	159.30
34170-121	Buffer pH 9.00	500ml	16.36	98.16
34170-258		4L	77.91	
34170-124	Buffer pH 10.00	500ml	18.67	112.02
34170-170		20L	171.91	
34170-133	Buffer pH 10.00 Blue	500ml	19.83	118.98
34170-161		20L	87.53	
34170-178		20 x 20ml	64.27	
34170-266	Buffer pH 11.00	500ml	22.95	
34170-268	Buffer pH 12.00	500ml	15.13	
34170-270		4L	72.18	

All VWR pH buffers are standardized to NIST Standard Reference Materials (SRM) and are accurate to +/-0.02 pH units. Our pH 4 and 7 buffers are accurate to +/-0.01 pH units. VWR buffers are available from pH 1.00 to 12.45 in colorless formulations. The pH 4, 7, and 10 buffers are available as either colorless or color coded in 500ml dispenser top poly bottles or in 4 or 20 liter Cubitainers. Color coded pH 4, 7, and 10 buffers are also available in 20ml poly lined single use sachets for lab or field use.

**VWR:**  
Providing value at competitive prices



**VWR**   
INTERNATIONAL

**vwr.com** 1.800.932.5000

Order from Over 750,000 Products

### Sales & Inventory Locations

#### Pacific Northwest

Anchorage, AK  
Boise, ID  
Brisbane, CA  
Seattle, WA  
Tualatin, OR

#### Southwest

Albuquerque, NM  
Aurora, CO  
Phoenix, AZ  
San Diego, CA  
San Dimas, CA

#### Midwest

Batavia, IL  
Indianapolis, IN  
Minneapolis, MN  
St. Louis, MO

#### Gulf

Houston, TX  
Lake Charles, LA  
Sugarland, TX

#### Northeast

Marlborough, MA  
Rochester, NY

#### Mid-Atlantic

Bridgeport, NJ  
Clarksburg, MD

#### Southeast

Miami, FL  
Morrisville, NC  
Oak Ridge, TN  
Suwanee, GA

#### Canada

Mississauga, Ontario  
Ville Mont-Royal,  
Québec

### Call Direct for Specialized Service Locations

#### International

911 Commerce Ct.  
Buffalo Grove, IL 60089  
847-520-6170  
fax: 847-229-0788

#### Puerto Rico

Carr. #869, Km. 1.5 M-4  
Royal Industrial Park  
Catano, PR 00962  
800-932-5000  
787-275-4500  
fax: 787-788-4320

#### Global and Corporate Accounts

1310 Goshen Pkwy.  
West Chester, PA 19380  
610-431-1700  
fax: 610-431-9174

#### Production Supplies and Services

975 Overland Ct.  
San Dimas, CA 91773  
877-VWR-PSCS  
(877-897-7727)  
fax: 888-890-9124

#### Furniture Division

3021 Gateway Drive  
Suite 280  
Irving, TX 75063  
888-624-2432  
fax: 972-753-1389

#### VWR Laboratory Distributors Services

911 Commerce Ct.  
Buffalo Grove, IL 60089  
800-873-8971  
fax: 847-229-0668

Prices are current at time of printing and are subject to change without notice. Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in U.S. dollars.

©2003 VWR International, Inc. All rights reserved. Printed in U.S.A.

09/03 (8m) Lit. No. 05321